

Anti- p24 antibody

Product Information

Catalog No.:	FNab09919
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5C6
IsoType:	IgG1
Storage:	0.01MPBS pH 7.4, -20°C for 12months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Recombinant Capsid protein p24
Synonyms:	None
Observed MW:	None
Uniprot ID :	P04591

Application

Reactivity:	HIV-1
Tested Application:	ELISA

FOR RESEARCH USE ONLY!

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com