

---

## anti- ZEB2 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab09954   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Immunogen information

|              |   |
|--------------|---|
| Immunogen:   | zinc finger E-box binding homeobox 2  |
| Synonyms:    | KIAA0569, SIP 1, SIP1, Smad interacting protein 1, SMADIP1, ZEB2, ZFHX1B, ZFX1B |
| Observed MW: | 136 kDa   |
| Uniprot ID : | O60315  |

### Application

|                       |   |
|-----------------------|---|
| Reactivity:           | Human, Mouse, Rat                                 |
| Tested Application:   | ELISA, WB, IHC, IF                                |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100 |
| Image:                |   |


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab09954( ZEB2 Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab09954( ZEB2 Antibody) at dilution of 1:1000