
anti- PPLA antibody

Product Information

Catalog No.:	FNab09994
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Cardiac phospholamban
Synonyms:	PLB, PLN, Cardiac phospholamban
Observed MW:	24 kDa
Uniprot ID :	P26678

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded rat heart tissue slide using FNab09994(CD31 Antibody) at dilution of 1:100



rat heart were subjected to SDS PAGE followed by western blot with FNab09993(CD31 Antibody) at dilution of 1:1000