

# anti- SULT1A1 antibody

#### **Product Information**

Catalog No.: FNab10011

Size: 100µg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. This gene encodes one of two phenol sulfotransferases with thermostable enzyme activity. Multiple alternatively spliced variants that encode two isoforms have been identified for this gene.

#### **Immunogen information**

Immunogen: sulfotransferase family, cytosolic, 1A, phenol-preferring, member 1

Synonyms: STP, STP1, ST1A1, Sulfotransferase 1A1

Observed MW: 33 kDa Uniprot ID: P50225

## **Application**

Reactivity: Mouse

Tested Application: ELISA, WB

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Recommended dilution: WB: 1:500 - 1:2000

Image:



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab10011(SULT1A1 antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>