

## anti- ELK3 antibody

### Product Information

Catalog No.:	FNab10053
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	ETS domain-containing protein Elk-3
Synonyms:	ETS-related protein ERP, Serum response factor accessory protein 2, SRF accessory protein 2, SAP-2, ETS-related protein NET
Observed MW:	44 kDa
Uniprot ID :	P41970

### Application

Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:500-1:2000; IHC: 1:50-1:500
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemical of paraffin-embedded human colon cancer using FNab10053(ELK3 antibody) at dilution of 1:100



K562 cells were subjected to SDS PAGE followed by western blot with FNab10053(ELK3 antibody) at dilution of 1:1000.