

# anti- OTX2 antibody

#### **Product Information**

Catalog No.: FNab10083

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

# **Background**

Probably plays a role in the development of the brain and the sense organs. Can bind to the BCD target sequence(BTS): 5'-TCTAATCCC-3'.

## **Immunogen information**

Immunogen: orthodenticle homeobox 2

Synonyms: Homeobox protein OTX2, MCOPS5, orthodenticle homeobox 2,

Orthodenticle homolog 2, OTX2

Observed MW: 35 kDa Uniprot ID: P32243

# Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:


1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded rat brain using FNab10083(OTX2 antibody) at dilution of 1:100



mouse liver were subjected to SDS PAGE followed by western blot with FNab10083(OTX2 antibody) at dilution of 1:1000