
anti- PGBM antibody

Product Information

Catalog No.:	FNab10133
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	heparan sulfate proteoglycan 2
Synonyms:	heparan sulfate proteoglycan 2, HSPG, HSPG2, Perlecan, PLC, PRCAN, SJA, SJS, SJS1
Observed MW:	469 kDa
Uniprot ID :	P98160

Application

Reactivity:	Human; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:100-1:400
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer using FNab10133(PGBM antibody) at dilution of 1:100