

anti- MKK7 antibody

Product Information

Catalog No.: FNab10158

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: mitogen-activated protein kinase kinase 7

Synonyms: JNK activating kinase 2, JNK kinase 2, JNKK 2, Jnkk2, MAP kinase

kinase 7, MAP2K7, MAPK/ERK kinase 7, MAPKK 7, MAPKK7, MEK

7, MEK7, MKK7, PRKMK7

Observed MW: 47 kDa, 51 kDa

Uniprot ID: O14733

Application

Reactivity: Human, rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded lung cancer using FNab10158(MKK7 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab10158(MKK7 antibody) at dilution of 1:1000