

anti- TEAD1 antibody

Product Information

Catalog No.:	FNab10169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Transcriptional enhancer factor TEF-1
Synonyms:	AA, NTEF 1, Protein GT IIC, REF1, TCF 13, TCF13, TEA domain family member 1, TEAD 1, TEAD1, TEF 1, TEF1, Transcription factor 13
Uniprot ID :	P28347

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10169 (TEAD1 antibody) at dilution of 1:200



Various lysates were subjected to SDS PAGE followed by western blot with FNab10169 (TEAD1 antibody) at dilution of 1:1000