

anti- MATK antibody

Product Information

Catalog No.: FNab10177

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: Megakaryocyte-associated tyrosine-protein kinase

Synonyms: CHK, CTK, HHYLTK, HYL, HYLTK, Lsk, MATK, Protein kinase HYL,

Tyrosine protein kinase CTK

Observed MW: 50-56 kDa Uniprot ID: P42679

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab10177(MATK antibody) at dilution of 1:200



K562 cells were subjected to SDS PAGE followed by western blot with FNab10177(MATK antibody) at dilution of 1:1000