

anti- SHBG antibody

Product Information

Catalog No.: FNab10201

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: Sex hormone-binding globulin

Synonyms: ABP, Androgen binding protein, SBP, sex hormone binding globulin,

Sex steroid binding protein, SHBG, TEBG

Observed MW: 32 kDa, 42 kDa, 44 kDa

Uniprot ID: P04278

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IHC, WB

Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human kidney cancer tissue slide using FNab10201(SHBG antibody) at dilution of 1:200



Rat testis were subjected to SDS PAGE followed by western blot with FNab10201(SHBG antibody) at dilution of 1:1000