

# anti- HMGB1 antibody

#### **Product Information**

Catalog No.: FNab10218

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen: High mobility group protein B1

Synonyms: high mobility group box 1, High mobility group protein 1, High mobility

group protein B1, HMG 1, HMG1, HMG3, HMGB1, SBP 1

Observed MW: 25 kDa Uniprot ID: P09429

## **Application**

Reactivity: Human

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:


#### 1

# Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Jurkat cells were subjected to SDS PAGE followed by western blot with FNab10218(HMGB1 antibody) at dilution of 1:1000