

# anti- CD280 antibody

### **Product Information**

Catalog No.: FNab10304

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Immunogen information**

Immunogen: MRC2

Synonyms: CD280, CLEC13E, ENDO180, UPARAP, MRC2

Observed MW: 190 kDa Uniprot ID: Q9UBG0

## **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500 - 1:2000

Image:


1

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab10304(CD280 antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>