

anti- 2019 nCOV N protein antibody

Product Information

Catalog No.: FNab10309

Size: 100μg Form: liquid

Purification: protein A purified

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 3C6
IsoType: IgG1

Storage: CB buffer, pH 7.5,-20°C for 12 months(Avoid repeated freeze / thaw

cycles.)

Immunogen information

Immunogen: 2019 nCOV N protein

Observed MW: 48 kDa

Application

Reactivity: 2019 nCOV

Tested Application: ELISA, WB, Colloidal Gold-Based

Recommended dilution: ELISA: Use at an assay dependent concentration (Suggestions: Use a

concentration of $0.1 - 1 \mu g/ml$)


Image:

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





2019 nCOV N protein were subjected to SDS PAGE followed by western blot with FNab10309 (anti- 2019 nCOV N protein Monoclonal antibody) at dilution of 1ug/ml