

anti- ARPC5L antibody

Product Information

Catalog No.:	FNab10385
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor(NPF) mediates the formation of branched actin networks.

Immunogen information

Immunogen:	actin related protein 2/3 complex, subunit 5-like
Synonyms:	ARC16 2, ARPC5L
Observed MW:	None
Uniprot ID :	Q9BPX5


Application

Reactivity:	Human, Rat		
Tested Application:	ELISA, IHC		
Recommended dilution: IHC: 1:50-1:200			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275





Immunohistochemistry of paraffin-embedded rat brain using FNab10385(ARPC5L antibody) at dilution of 1:100

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com