

anti- HSP27 antibody

Product Information

Catalog No.:	FNab10398	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

HSPB1, also known as heat shock protein 27(Hsp27), belongs to the small heat shock protein family which is induced in response to environmental challenges or/and developmental transitions. It is also an anti-apoptotic protein that plays crucial roles in tumorigenesis and cell survival and is reported to be an independent prognosis marker for cancer. Recently HSPB1 has been found to be a valuable marker for melanoma. In addition to the predicted 27 kDa, an extra 50-55 kDa representing dimeric form of HSPB1 may also be observed.(21353161)

Immunogen information

Immunogen:	heat shock 27kDa protein 1
Synonyms:	Heat shock protein beta 1, HSP 27, Hsp25, HSP27, HSP28, HSPB1
Observed MW:	19-27 kDa
Uniprot ID :	P04792

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275

www.fn-test.com



Recommended dilution:WB: 1:500-1:3000 Image:



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04044(HSPB1 antibody) at dilution of 1:1000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com