

## anti- 2019 nCOV S1 protein antibody

## **Product Information**

Catalog No.:	FNab10400
Size:	100µg
Form:	liquid
Purification:	protein A purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Human
Clonality:	monoclonal
Clone ID:	7H4
IsoType:	IgG
Storage:	PBS pH 7.4,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Immunogen information	
Immunogen:	S1-RBD

## Application

Reactivity:	2019 nCOV S protein
Tested Application:	ELISA, WB
Recommended dilution: ELISA: Use at an assay dependent concentration (Suggestions: Use concentration of 0.1 - 1 µg/ml);	
	WB: 1:1000-1:10000

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

1