

anti- ADIPOR1 antibody

Product Information

Catalog No.:	FNab10402	
Size:	100µg	
Form:	liquid	
Purification:	Protein A+G purification	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Mouse	
Clonality:	monoclonal	
Clone ID:	3C3	
IsoType:	IgG1	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Immunogen information

Immunogen:	adiponectin receptor 1
Synonyms:	ACDCR1, Adiponectin receptor 1, Adiponectin receptor protein 1, ADIPOR1, CGI 45, CGI45, PAQR1, TESBP1A
Observed MW:	43 kDa
Uniprot ID :	Q96A54


Application

Reactivity:	Human		
Tested Application:	ELISA, WB		
Recommended dilution: WB: 1:500 - 1:3000			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com