

anti- GFPT1 antibody

Product Information

Catalog No.: FNab10406

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

This gene encodes the first and rate-limiting enzyme of the hexosamine pathway and controls the flux of glucose into the hexosamine pathway. The product of this gene catalyzes the formation of glucosamine 6-phosphate.

Immunogen information

Immunogen: glutamine-fructose-6-phosphate transaminase 1

Synonyms: GFAT, GFPT

Observed MW:

Uniprot ID: Q06210

Application

Reactivity: Human

Tested Application: ELISA, IHC

Recommended dilution: IHC: 1:50 - 1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human liver cancer using FNab10406(GFPT1 antibody) at dilution of 1:100

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>