
anti- Caveolin-1 antibody

Product Information

Catalog No.:	FNab10420
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	caveolin 1, caveolae protein, 22kDa
Synonyms:	BSCL3, CAV, CAV1, Caveolin 1, Caveolin-1, CGL3, MSTP085, VIP21
Observed MW:	50 kDa
Uniprot ID :	Q03135

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab10420(Caveolin-1 Antibody) at dilution of 1:200



Various lysates were subjected to SDS PAGE followed by western blot with FNab10420(CAV1 antibody) at dilution of 1:1000