

anti- GATA4 antibody

Product Information

Catalog No.: FNab10433

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: GATA binding protein 4

Synonyms: GATA binding factor 4, GATA binding protein 4, GATA4,

Transcription factor GATA 4

Observed MW: 45-55kd Uniprot ID: P43694

Application

Reactivity: Human, Rat Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Various lysates were subjected to SDS PAGE followed by western blot with FNab10433(GATA4 antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>