

anti- A β 42 antibody

Product Information

Catalog No.:	FNab10440
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Amyloid-beta protein 42
Synonyms:	Amyloid-beta protein 42, Beta-APP42, Abeta42
Observed MW:	72-86kd, 100-140 kDa
Uniprot ID :	P05067

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain tissue slide using FNab10440(A β 42 Antibody) at dilution of 1:100



293T cells were subjected to SDS PAGE followed by western blot with FNab10440(A β 42 Antibody) at dilution of 1:1000.