

# anti- BIRC8 antibody

## **Product Information**

Catalog No.:	FNab10464
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

## Background

Protects against apoptosis mediated by BAX.

# **Immunogen information**

Immunogen:	baculoviral IAP repeat-containing 8
Synonyms:	BIRC8, hILP2, IAP like protein 2, ILP 2, ILP2
Observed MW:	30-40 kDa
Uniprot ID :	Q96P09

# Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, IHC, WB		
Recommended dilution: WB: 1:200-1:2000; IHC: 1:50-1:200			
Image:			

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :( 0086)027-87384275


www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer using FNab10464(BIRC8 antibody) at dilution of 1:100

various lysates were subjected to SDS PAGE followed by western blot with FNab10464(BIRC8 antibody) at dilution of 1:1000



2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

www.fn-test.com