

anti- GSTO1 antibody

Product Information

Catalog No.: FNab10470

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Exhibits glutathione-dependent thiol transferase and dehydroascorbate reductase activities. Has S-(phenacyl)glutathione reductase activity. Has also glutathione S-transferase activity. Participates in the biotransformation of inorganic arsenic and reduces monomethylarsonic acid(MMA) and dimethylarsonic acid.

Immunogen information

Immunogen: glutathione S-transferase omega 1

Synonyms: DKFZp686H13163, GSTO 1, GSTO1, GSTTLp28, P28

Observed MW: 28kd Uniprot ID: P78417

Application

Reactivity: Human, Rat
Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Image:



rat brain were subjected to SDS PAGE followed by western blot with FNab10470(GSTO1 antibody) at dilution of 1:500

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>