

anti-NCR1 antibody

Product Information

Catalog No.:	FNab10472
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

The natural cytotoxicity receptors (NCRs) are a recently characterized family of Ig-like activation receptors that appear to be major triggering receptors in tumor cell recognition. NCR1 is a glycoprotein that has two extracellular Ig-like domains followed by a ~40 amino acid residue stalk region, a type I transmembrane domain, and a short cytoplasmic tail. NCR1 has been shown to represent a novel NK cell-specific molecule involved in human NK cell activation. NCR1 has been implicated in NK cell-mediated lysis of several autologous tumor cells and pathogen-infected cell lines.

Immunogen information

Immunogen:	Natural cytotoxicity triggering receptor 1
Synonyms:	CD335, LY94, NK-p46, NKP46, NCR1
Observed MW:	46 kDa
Uniprot ID :	O76036

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:50-1:200

Image:



Immunohistochemistry of paraffin-embedded human skin cancer slide using FNab10472(NCR1 Antibody) at dilution of 1:200



A431 cells were subjected to SDS PAGE followed by western blot with FNab10472(NCR1 Antibody) at dilution of 1:1000