

anti- CPT1C antibody

Product Information

Catalog No.:	FNab10524
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play a role in lipid metabolic process.

Immunogen information

Immunogen:	carnitine palmitoyltransferase 1C
Synonyms:	CATL1, CPT 1C, CPT1 B, CPT1C, CPT1P, CPTI B, CPTIC
Uniprot ID :	Q8TCG5

Application

Reactivity:	Human, Rat
Tested Application:	IHC
Recommended dilution:IHC:	1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain using FNab10524(CPT1C antibody) at dilution of 1:50