

## anti- SCN7A antibody

### Product Information

Catalog No.:	FNab10527
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Peptide
Synonyms:	SCN6A, SCN7A, sodium channel, voltage-gated, type VII, alpha
Observed MW:	193kd
Uniprot ID :	Q01118

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded mouse heart tissue slide using FNab10527(SCN7A Antibody) at dilution of 1:100



mouse heart were subjected to SDS PAGE followed by western blot with FNab10527(SCN7A antibody) at dilution of 1:1000