

anti- Phospho- β -Catenin(S33/S37/T41) Antibody

Product Information

Catalog No.:	FNab10701
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A synthetic phosphorylated peptide around S33 & T41 of human beta Catenin
Synonyms:	Phospho- β -Catenin (Ser33/37/Thr41), Phospho- β -Catenin-S33/S37/T41
Observed MW:	90 kDa
Uniprot ID :	P35222

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



293T cells were subjected to SDS PAGE followed by western blot with FNab10701 at dilution of 1:1000