

anti- LGP2 antibody

Product Information

Catalog No.:	FNab10709
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

DHX58(Probable ATP-dependent RNA helicase DHX58) is also named as D11LGP2E, LGP2 and belongs to the RLR subfamily.It acts as a regulator of DDX58/RIG-I and IFIH1/MDA5 mediated antiviral signaling.

Immunogen information

Immunogen:	DEXH(Asp-Glu-X-His) box polypeptide 58
Synonyms:	D11LGP2, D11lgp2e, DHX58, LGP2, Protein D11Lgp2 homolog, RIG I like receptor LGP2, RLR
Observed MW:	70-75 kDa, 35kDa
Uniprot ID :	Q96C10

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10709(LGP2 antibody) at dilution of 1:100



Mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab10709(LGP2 antibody) at dilution of 1:1000