

anti- MICB antibody

Product Information

Catalog No.:	FNab10743	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

Seems to have no role in antigen presentation. Acts as a stress-induced self-antigen that is recognized by gamma delta T cells. Ligand for the KLRK1/NKG2D receptor. Binding to KLRK1 leads to cell lysis.

Immunogen information

Immunogen:	MHC class I polypeptide-related sequence B	
Synonyms:	MICB, PERB11.2, MHC class I polypeptide-related sequence B	
Observed MW:	40-60 kDa	
Uniprot ID :	Q29980	

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, IHC, WB		
Recommended dilution:WB: 1:500-1:2000; IHC: 1:20-1:200			
Image:			

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206) Tel :(0086)027-87384275

www.fn-test.com





Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab10743(MICB antibody) at dilution of 1:100



L6 cell lysates were subjected to SDS PAGE followed by western blot with FNab10743(MICB antibody) at dilution of 1:1000

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

www.fn-test.com