

## anti- CD3 gamma antibody

### Product Information

Catalog No.:	FNab10749
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	CD3g molecule, gamma(CD3-TCR complex)
Synonyms:	CD3 gamma, CD3G, CD3-gamma, T cell receptor T3 gamma chain, T3G
Observed MW:	20-25 kDa
Uniprot ID :	P09693

### Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Raji cell lysates were subjected to SDS PAGE followed by western blot with FNab10749(CD3E antibody) at dilution of 1:1000