

anti- VNN1 antibody

Product Information

Catalog No.:	FNab10796
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Amidohydrolase that hydrolyzes specifically one of the carboamide linkages in D-pantetheine thus recycling pantothenic acid(vitamin B5) and releasing cysteamine.

Immunogen information

Immunogen:	vanin 1
Synonyms:	HDLCQ8, Pantetheinase, Pantetheine hydrolase, Tiff66, vanin 1, VNN1
Observed MW:	55-72 kDa
Uniprot ID :	O95497

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



A431 cells were subjected to SDS PAGE followed by western blot with FNab10796(VNN1 antibody) at dilution of 1:1000