

---

## anti- CD4 antibody

### Product Information

Catalog No.:	FNab10808
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	5H1
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	CD4 molecule
Synonyms:	CD4, CD4 molecule, CD4mut
Uniprot ID :	P01730

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution: IHC:	1:200 - 1:1000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human liver using FNab10808(CD4 antibody) at dilution of 1:500