
anti- progesterone antibody

Product Information

Catalog No.:	FNab10809
Size:	100µg
Form:	liquid
Purification:	protein A purified
Purity:	≥95% as determined by SDS-PAGE
Host:	mouse
Clonality:	monoclonal
Clone ID:	1G8
IsoType:	IgG1
Storage:	PBS pH 7.4 , 2 - 8° C for short time, -20° C for 12 months. (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Progesterone-11-BSA
Synonyms:	Progesterone, prog, PRG, PGN

Application

Reactivity:	universal
Tested Application:	ELISA, Lateral Flow, CLIA

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com