

## anti- FLNA antibody

### Product Information

Catalog No.:	FNab10810
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Filamin-A
Synonyms:	ABP 280, ABPX, Actin binding protein 280, Alpha filamin, DKFZp434P031, Filamin 1, Filamin A, FLN, FLN A, FLN1, FLNA, FMD, MNS, NHBP, Non muscle filamin, OPD, OPD1, OPD2
Observed MW:	300 kDa
Uniprot ID :	P21333

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast cancer using FNab10810(FLNA antibody) at dilution of 1:100



Hela cell lysates were subjected to SDS PAGE followed by western blot with FNab10810(FLNA antibody) at dilution of 1:1000